

M400 2-wire Transmitter Series for pH/ORP, Oxygen, Conductivity and Dissolved Carbon Dioxide

Versatile, intelligent transmitters for harsh conditions

Technical Data



Short description

The M400 2-wire transmitter series is the state-of-the-art transmitter for most demanding conditions in hazardous and non-hazardous area applications. The transmitter series features advanced ISM technology and covers pH/ORP, Oxygen, Dissolved Carbon Dioxide and Conductivity measurements. Thanks to the mixed-mode input functionality, the M400 accepts any analog or ISM sensor of your choice. The M400 is a single-channel, multi-parameter unit. The same unit can handle different parameters such as pH/ORP, Oxygen (for measurement of dissolved oxygen or in gas), Ozone, Dissolved Carbon Dioxide or Conductivity, depending on the type you choose.

Features

- IECEx/ATEX/FM/NEPSI approved version
- Advanced ISM functionalities
- Mixed-mode input (analog or ISM sensors accepted)
- Communication protocols: 4 to 20 mA (with HART®), FOUNDATION fieldbus®, PROFIBUS PA®
- Multi-parameter unit
- Compatible with optical dissolved oxygen sensors
- Dynamic Lifetime Indication
- Adaptive Calibration Timer
- Time to Maintenance
- Analog 4 to 20 mA input signal (for pressure compensation)
- IP66/NEMA 4X rated
- PID controller
- Quick setup mode
- 8 languages: English, German, French, Italian, Spanish, Portuguese, Russian and Japanese

ISM®

HART
COMMUNICATION PROTOCOL

Fieldbus
Foundation

PROFIBUS




Contents

Dimensions and installation drawings	2
Specifications	3
Terminal block definitions	10
Ordering information	14


METTLER TOLEDO

Assembly



- 1 5 pieces M20 × 1.5 cable glands
- 2 2 pieces plastic plugs
- 3 3 pieces screws

Dimension drawings



pH/ORP (incl. pH/pNa)

Measurement parameters	pH, mV and temperature
pH display range	-2.00 to +20.00 pH
pH resolution	Auto/0.001/0.01/0.1/1 (can be selected)
pH accuracy ¹⁾	Analog: ±0.02 pH
mV range	-1500 to +1500 mV
mV resolution	Auto/0.001/0.01/0.1/1 mV (can be selected)
mV accuracy ¹⁾	Analog: ±1 mV
Temperature input ²⁾	Pt1000/Pt100/NTC30K
Temperature measuring range	-30 to 130 °C (-22 to 266 °F)
Temperature resolution	Auto/0.001/0.01/0.1/1 (can be selected)
Temperature accuracy ¹⁾	Analog: ±0.25 K in the range of -10 to +150 °C (±0.45 °F in the range of +14 to +302 °F)
Temperature repeatability ¹⁾	±0.13 K (±0.23 °F)
Temperature compensation	Automatic/Manual
Max. sensor cable length	<ul style="list-style-type: none"> • Analog: 10 to 20 m (33 to 65 ft) depending on sensor • ISM: 80 m (260 ft)
Calibration	1-point (offset), 2-point (slope and offset) or Process (offset)

1) ISM input signal causes no additional error.
2) Not required on ISM sensors

Amperometric oxygen

Measurement parameters	<ul style="list-style-type: none"> • Dissolved oxygen: Saturation or concentration and temperature • Oxygen in gas: Concentration and temperature
Current range	Analog: 0 to -7000 nA
Oxygen measuring ranges, dissolved oxygen	<ul style="list-style-type: none"> • Saturation: 0 to 500% air, 0 to 200% O₂ • Concentration: 0 ppb (µg/L) to 50.00 ppm (mg/L)
Oxygen measuring ranges, oxygen in gas	0 to 9999 ppm O ₂ gas, 0 to 100 vol % O ₂
Oxygen accuracy, dissolved oxygen ¹⁾	<ul style="list-style-type: none"> • Saturation: ±0.5% of the measured value or ±0.5%, depending on which is larger • Concentration at high values: ±0.5% of the measured value or ±0.050 ppm/±0.050 mg/L, depending on which is larger • Concentration at low values: ±0.5% of the measured value or ±0.001 ppm/±0.001 mg/L, depending on which is larger • Concentration at traces values: ±0.5% of the measured value or ±0.100 ppb/±0.1 µg/L, depending on which is larger
Oxygen accuracy, oxygen in gas ¹⁾	<ul style="list-style-type: none"> • ±0.5% of the measured value or ±5 ppb, depending on which is larger for ppm O₂ gas • ±0.5% of the measured value or ±0.01%, depending on which is larger for vol % O₂
Resolution current ¹⁾	Analog: 6 pA
Polarization voltage	<ul style="list-style-type: none"> • Analog: -1000 to 0 mV (configurable) • ISM: -550 mV or -674 mV (configurable)
Temperature input	NTC 22 kΩ, Pt1000, Pt100
Temperature compensation	Automatic
Temperature measuring range	-10 to +80 °C (+14 to +176 °F)
Temperature accuracy	±0.25 K in the range of -10 to +80 °C (+14 to +176 °F)
Max. sensor cable length	<ul style="list-style-type: none"> • Analog: 20 m (65 ft) • ISM: 80 m (260 ft)
Calibration	1-point (slope and offset) or Process (slope and offset)

1) ISM input signal causes no additional error.

Optical dissolved oxygen

Measurement parameters	Dissolved oxygen (DO) saturation or concentration and temperature
DO concentration range	0.1 ppb (µg/L) to 50.00 ppm (mg/L)
DO saturation range	0 to 500% air, 0 to 100% O ₂
DO resolution	Auto/0.001/0.01/0.1/1 (can be selected)
DO accuracy	±1 digit
Temperature measuring range	-30 to +150 °C (-22 to +302 °F)
Temperature resolution	Auto/0.001/0.01/0.1/1 (can be selected)
Temperature accuracy	±1 digit
Temperature repeatability	±1 digit
Temperature compensation	Automatic
Max. sensor cable length	15 m (50 ft)
Calibration	1-point (depending on sensor model), 2-point or Process

Dissolved carbon dioxide (CO₂ low)

Measurement parameters	Dissolved carbon dioxide and temperature
CO ₂ measuring ranges	<ul style="list-style-type: none"> • 0 to 5000 mg/L • 0 to 200 %sat • 0 to 1500 mm Hg • 0 to 2000 mbar • 0 to 2000 hPa
CO ₂ accuracy	±1 digit
CO ₂ resolution	Auto/0.001/0.01/0.1/1 (can be selected)
mV range	-1500 to +1500 mV
mV resolution	Auto/0.01/0.1/1 mV
mV accuracy	±1 digit
Total pressure range (TotPres)	0 to 4000 mbar
Temperature input	Pt1000/NTC22K
Temperature measuring range	0 to +60 °C (+32 to +140 °F)
Temperature resolution	Auto/0.001/0.01/0.1/1, (can be selected)
Temperature accuracy	±1 digit
Temperature repeatability	±1 digit
Max. sensor cable length	80 m (260 ft)
Calibration	1-point (offset), 2-point (slope and offset) or Process (offset)

Thermal conductivity CO₂ hi (InPro 5500i) (PROFIBUS PA® only)

CO ₂ measuring ranges	<ul style="list-style-type: none"> • 0 to 10 bar p (CO₂)/0 to 145 psi p (CO₂) • 0 to 15 g/L • 0 to 7 V/V CO₂
Accuracy in fluids ¹⁾	<ul style="list-style-type: none"> • ±1 % of reading (within ±5 % of calibration temperature) • ±2 % of reading over temperature range 0 to +50 °C (+32 to +122 °F)

1) Complete loop of sensor and transmitter

Dissolved ozone (PROFIBUS PA® only)

Measurement parameters	Concentration and temperature
Display range for current	0 to –900 nA
Ozone measuring range	Concentration 0.1 ppb (µg/L) to 5.00 ppm (mg/L) O ₃
Ozone accuracy	± 1 digit
Resolution current	± 1 digit
Temperature compensation	Automatic
Temperature measuring range	0 to + 50 °C (+32 to + 122 °F)
Temperature resolution	Auto/0.001/0.01/0.1/1 (can be selected)
Temperature accuracy	± 1 digit
Max. sensor cable length	80 m
Calibration	1-point ZeroPt or Process (ZeroPt and slope)

Conductivity 2-e/4-e

Measurement parameters	Conductivity/resistivity and temperature
Conductivity ranges	0.02 to 2,000 µS/cm (500 Ω × cm to 50 MΩ × cm)
2-electrode sensor	C = 0.01 0.002 to 200 µS/cm (5000 Ω × cm to 500 MΩ × cm)
	C = 0.1 0.02 to 2000 µS/cm (500 Ω × cm to 50 MΩ × cm)
	C = 1 15 to 4000 µS/cm
	C = 3 15 to 12,000 µS/cm
	C = 10 10 to 40,000 µS/cm (25 Ω × cm to 100 kΩ × cm)
Conductivity ranges	0.01 to 650 mS/cm (1.54 Ω × cm to 0.1 MΩ × cm)
4-electrode sensor	
Display range for 2-e sensor	0 to 40,000 mS/cm (25 Ω × cm to 100 MΩ × cm)
Display range for 4-e sensor	0.01 to 650 mS/cm (1.54 Ω × cm to 0.1 MΩ × cm)
Chemical concentration curves	NaCl: 0–26% @ 0 °C to 0–28% @ +100 °C NaOH: 0–12% @ 0 °C to 0–16% @ + 40 °C to 0–6% @ +100 °C HCl: 0–18% @ –20 °C to 0–18% @ 0 °C to 0–5% @ +50 °C HNO ₃ : 0–30% @ –20 °C to 0–30% @ 0 °C to 0–8% @ +50 °C H ₂ SO ₄ : 0–26% @ –12 °C to 0–26% @ + 5 °C to 0–9% @ +100 °C H ₃ PO ₄ : 0–35% @ + 5 °C to + 80 °C User-defined concentration table (5 × 5 matrix)
TDS ranges	NaCl, CaCO ₃
Cond/Res accuracy ¹⁾	Analog: ±0.5 % of reading or 0.25 Ω, whichever is greater, up to 10 MΩ-cm
Cond/Res repeatability ¹⁾	Analog: ±0.25% of reading or 0.25 Ω, whichever is greater
Cond/Res resolution	Auto/0.001/0.01/0.1/1 (can be selected)
Temperature input	Pt1000/Pt100/NTC22K
Temperature measuring range	–40 to + 200 °C (–40 to + 392 °F)
Temperature resolution	Auto/0.001/0.01/0.1/1 (can be selected)
Temperature accuracy	• ISM: ± 1 digit • Analog: ±0.25 K (±0.45 °F) within –30 to + 150 °C (–22 to + 302 °F); ±0.50 K (±0.90 °F) outside
Temperature repeatability ¹⁾	±0.13 K (±0.23 °F)
Max. sensor cable length	• ISM: 80 m (260 ft) • Analog: 61 m (200 ft); with 4-e sensors: 15 m (50 ft)
Calibration	1-point, 2-point or Process

1) ISM input signal causes no additional error.

Specifications

M400 2-wire transmitter series

Inductive Conductivity (M400 Cond Ind transmitter only)

Measurement parameters	Conductivity and temperature
Display range	0 to 2,000 mS/cm
Chemical concentration curves	NaCl: 0–26 % @ 0 °C to 0 – 28 % @ +100 °C NaOH-1: 0–13 % @ 0 °C to 0 – 24 % @ +100 °C NaOH-3: 15–50 % @ 0 °C to 35 – 50 % @ +100 °C HCl-1: 0–18 % @ –20 °C to +50 °C HCl-2: 22–39 % @ –20 °C to +50 °C HNO ₃ -1: 0–30 % @ –20 °C to +50 °C HNO ₃ -2: 35–96 % @ –20 °C to +50 °C H ₂ SO ₄ -1: 0–26 % @ –12 °C to 0–37 % @ +100 °C H ₂ SO ₄ -2: 28–88 % @ 0 °C to 39–88 % @ +95 °C H ₂ SO ₄ -3: 94–99 % @ –12 °C to 89–99 % @ +95 °C H ₃ PO ₄ : 0–35 % @ + 5 °C to + 80 °C User-defined concentration table (5 × 5 matrix)
TDS ranges	NaCl, CaCO ₃
Conductivity accuracy	± 1.0 % of reading or ± 0.005 mS/cm
Conductivity repeatability	± 1.0 % of reading or ± 0.005 mS/cm
Conductivity resolution	Auto/0.001/0.01/0.1/1 (can be selected)
Temperature input	Pt1000/Pt100/NTC22K
Temperature measuring range	–40 to +200 °C (–40 to +392 °F)
Temperature resolution	Auto/0.001/0.01/0.1/1 (can be selected)
Temperature accuracy	± 0.25 K (± 0.45 °F) within –30 to +150 °C (–22 to +302 °F); ± 0.50 K (± 0.90 °F) outside
Temperature repeatability	± 0.13 K (± 0.23 °F)
Max. sensor cable length	10 m (32.8 ft)
Calibration	1-point, Zero point or Process

General electrical specifications

Display	Backlit LCD, 4 lines
Running capacity	Ca. 4 days
Keypad	5 tactile feedback keys
Languages	8 (English, German, French, Italian, Spanish, Portuguese, Russian and Japanese)
Connection terminals	Spring cage terminals, appropriate for wire cross section 0.2 to 1.5 mm ² (AWG 16 – 24)
Analog input	4 to 20 mA (for pressure compensation)

Specification for 4 to 20 mA (with HART ®)

Supply voltage	14 to 30 V DC
Number of outputs (analog)	2
Current outputs	Loop current 4 ... 20 mA, galvanically isolated up to 60 V from input and from earth/ground, protected against wrong polarity, feeding voltage 14 to 30 V DC
Measurement error through analog outputs	<±0.05 mA over 1 to 20 mA range
Analog output configuration	Linear
PID process controller	Pulse length, pulse frequency
Hold input/Alarm contact	Yes/Yes (alarm delay 0 to 999 s)
Digital outputs	2 open collector (OC), 30 V DC, 100 mA, 0.9 W
Digital input	<ul style="list-style-type: none"> • M400/2H, M400/2XH, M400G/2H, M400G/2XH: 2 • M400/2XH Cond Ind: 1 • Galvanically isolated up to 60 V from output, analog input and ground/earth with switching limits 0.00 V DC to 1.00 V DC inactive 2.30 V DC to 30.00 V DC active
Alarm output delay	0 to 999 s

Specification for FOUNDATION fieldbus®

Supply voltage	<ul style="list-style-type: none"> • Non hazardous area (Non-IS): 9 to 32 V DC • Linear Barrier: 9 to 24 V DC • FISCO: 9 to 17.5 V DC
Current	22 mA
Max. current in case of fault (FDE)	<28 mA
Number of current inputs	1 for pressure compensation
Physical interface	According to IEC 61158-2
Transfer rate	31.25 kbit/s
Profile	FF_H1 (Foundation fieldbus)
Communication protocol	FF-816
ITK version	6.0.1
Manufacturer ID (DEV_TYPE)	0x465255
FF Type (DEV_REV)	1
FF communication model	<ul style="list-style-type: none"> • 1 Resource Block • 1 Physical Block • 2 Transducer Blocks (General and Sensor) • 4 Analog Input Blocks • 1 Analog Output Block • 2 Discrete Input Blocks • 2 Discrete Output Blocks

PROFIBUS PA® specifications

Supply voltage	Non hazardous area (Non-IS): 9 to 32 V DC Linear Barrier: 9 to 24 V DC FISCO: 9 to 17.5 V DC
Current consumption	22 mA
Current consumption on error	< 28 mA
Number of current inputs	1 for pressure compensation
Profile	PROFIBUS PA 3.02
PA communication model	<ul style="list-style-type: none"> • 1 Resource Block • 1 Physical Block • 1 Analyser Transducer Block (Sensor Block) • 4 Analog Input Blocks • 1 Analog Output Block • 2 Discrete Input Blocks • 2 Discrete Output Blocks

Environmental specifications

Storage temperature	-40 to +70 °C (-40 to +158 °F)
Ambient temperature operating range	-20 to +60 °C (-4 to +140 °F)
Relative humidity	0 to 95 % non-condensing
EMC	According to EN 61326-1 (general requirements) Emission: Class B, Immunity: Class A
Certificates and approvals	<p>M400/2H</p> <ul style="list-style-type: none"> • cFMus Class I, Division 2, Groups A, B, C, D T4A • cFMus Class I, Zone 2, Groups IIC T4 <hr/> <p>M400/2XH, M400G/2XH, M400/2XH Cond Ind</p> <ul style="list-style-type: none"> • ATEX/IECEX Zone 1 Ex ib [ia Ga] IIC T4 Gb • ATEX/IECEX Zone 21 Ex ib [ia Da] IIIC T80°C Db IP66 • cFMus Class I, Division 1, Groups A, B, C, D T4 • cFMus Class II, Division 1, Groups E, F, G • cFMus Class III • cFMus Class I, Zone 0, AEx ia IIC T4 Ga • NEPSI EX Zone <hr/> <p>M400FF</p> <ul style="list-style-type: none"> • ATEX/IECEX Zone 1 Ex ib [ia Ga] IIC T4 Gb • cFMus Class I, Division 1, Groups A, B, C, D T4A • NEPSI EX Zone <hr/> <p>M400PA</p> <ul style="list-style-type: none"> • ATEX/IECEX Zone 1 Ex ib [ia Ga] IIC T4 Gb • cFMus Class I, Division 1, Groups A, B, C, D T4A • NEPSI EX Zone
CE mark	The measuring system is in conformity with the statutory requirements of the EC Directives. METTLER TOLEDO confirms successful testing of the device by affixing to it the CE mark.


Specifications

M400 2-wire transmitter series

Mechanical specifications

Dimensions	Housing –	144 × 144 × 116 mm
	Height × Width × Depth	(5.7 × 5.7 × 4.6 inch)
	Front bezel –	150 × 150 mm
	Height × Width	(5.9 × 5.9 inch)
	Max. depth – panel mounted	87 mm (excludes plug-in connectors)
Weight		1.50 kg (3.3 lb)
Material		Aluminum die cast
Enclosure rating		IP 66/NEMA4X

Terminal block (TB) definitions



- 1 TB1 – Input and output analog signal
- 2 TB2 – Sensor signal

TB1 terminal definition

4 to 20 mA (with HART®)

Terminal	Description
1	DI1+
2	DI1-
3	DI2+
4	DI2-
5	Not used
6	OC1+
7	OC1-
8	OC2+
9	OC2-
10	AO1+/HART
11	AO1-/HART
12	AO2+
13	AO2-
14	not used
15	⊥

Foundation fieldbus®

Terminal	Description
1	Not available
2	Not available
3	Not available
4	Not available
5	Not available
6	Not available
7	Not available
8	Not available
9	Not available
10	+FF-H1
11	-FF-H1
12	+FF-H1
13	-FF-H1
14	not used
15	⊥

PROFIBUS PA®

Terminal	Description
1	Not available
2	Not available
3	Not available
4	Not available
5	Not available
6	Not available
7	Not available
8	Not available
9	Not available
10	+PA
11	-PA
12	+PA
13	-PA
14	not used
15	⊥

TB2 terminal definition – Analog sensors (except M400 Cond Ind transmitter)**Conductivity 2-e/4-e**

Terminal	Function	Color
A	Cnd inner1 ¹⁾	White
B	Cnd outer1 ¹⁾	White/blue
C	Cnd outer1	–
D	Not used	–
E	Cnd outer2	–
F	Cnd inner2 ²⁾	Blue
G	Cnd outer2 (GND) ²⁾	Black
H	Not used	–
I	RTD ref/GND	Bare shield
J	RTD sense	Red
K	RTD	Green
L	Not used	–
M	Not used	–
N	Not used	–
O	Not used	–
P	Not used	–
Q	Not used	–

1) For third party Conductivity 2-e sensors a jumper between A and B may be required.

2) For third party Conductivity 2-e sensors a jumper between F and G may be required.

pH/ORP

Terminal	pH		Redox (ORP)	
	Function	Color ¹⁾	Function	Color
A	Glass	Transparent	Platinum	Transparent
B	Not used	–	–	–
C	Not used	–	–	–
D	Not used	–	–	–
E	Reference	Red	Reference	Red
F	Reference ²⁾	–	Reference ²⁾	–
G	Solution GND ²⁾	Blue ³⁾	Solution GND ²⁾	–
H	Not used	–	–	–
I	RTD ref/GND	White	–	–
J	RTD sense	–	–	–
K	RTD	Green	–	–
L	Not used	–	–	–
M	Shield (GND)	Green/yellow	Shield (GND)	Green/yellow
N	Not used	–	–	–
O	Not used	–	–	–
P	Not used	–	–	–
Q	Not used	–	–	–

1) Grey wire not used.

2) Install jumper between F and G for ORP sensors and pH electrodes without SG.

3) Blue wire for electrode with SG.

TB2 terminal definition – Analog sensors (continuation) (except M400 Cond Ind transmitter)**Amperometric oxygen**

Terminal	Function	InPro 6800(G)	InPro 6900	InPro 6950
		Color	Color	Color
A	Not used	–	–	–
B	Anode	Red	Red	Red
C	Anode	– 1)	– 1)	–
D	Reference	– 1)	– 1)	Blue
E	Not used	–	–	–
F	Not used	–	–	–
G	Guard	–	Grey	Grey
H	Cathode	Transparent	Transparent	Transparent
I	NTC ref (GND)	White	White	White
J	Not used	–	–	–
K	NTC	Green	Green	Green
L	Not used	–	–	–
M	Shield (GND)	Green/yellow	Green/yellow	Green/yellow
N	Not used	–	–	–
O	Not used	–	–	–
P	+Ain ²⁾	–	–	–
Q	–Ain ²⁾	–	–	–

1) Install jumper between C and D for InPro 6800(G) and InPro 6900.

2) 4 to 20 mA signal for pressure compensation

TB2 terminal definition – Analog sensors (M400 Cond Ind transmitter only)**Inductive conductivity**

Terminal	Function	Color	
		InPro 7259 ST, InPro 7250 PFA	InPro 7250 HT
A	Not used	–	–
B	Not used	–	–
C	Not used	–	–
D	Send High	Blue	Black or Transparent
E	Send Low	Brown	Violet
F	Shield (GND)	Green-Yellow	Green-Yellow
G	Receive Low	Red	Yellow
H	Receive High	Black or Transparent	Red
I	RTD	White	White
J	RTD sense	Grey	Grey
K	RTD	Green	Green
L – Q	Not used	–	–

TB2 terminal definition – ISM sensors (except M400 Cond Ind transmitter)**pH, Amperometric oxygen, Ozone¹⁾, Conductivity 4-e, Dissolved CO₂ low**

Terminal	Function	Color
A	Not used	–
B	Not used	–
C	Not used	–
D	Not used	–
E	Not used	–
F	Not used	–
G	Not used	–
H	Not used	–
I	Not used	–
J	Not used	–
K	Not used	–
L	1-wire	Transparent (cable core)
M	GND	Red (shield)
N	RS485-B	–
O	RS485-A	–
P	+Ain ²⁾	–
Q	–Ain ²⁾	–

1) PROFIBUS PA® only

2) Only for Oxygen sensors: 4 to 20 mA signal for pressure compensation

Optical dissolved oxygen, CO₂ hi (InPro 5500 i) ¹⁾

Terminal	Optical dissolved oxygen with VP8 cable ²⁾		Optical dissolved oxygen with other cables ³⁾	
	Function	Color	Function	Color
A	Not used	–	Not used	–
B	Not used	–	Not used	–
C	Not used	–	Not used	–
D	Not used	–	Not used	–
E	Not used	–	Not used	–
F	Not used	–	Not used	–
G	Not used	–	Not used	–
H	Not used	–	Not used	–
I	Not used	–	D_GND (shield)	Yellow
J	Not used	–	Not used	–
K	Not used	–	Not used	–
L	Not used	–	Not used	–
M	D_GND (shield)	Green / yellow	D_GND (shield)	Grey
N	RS485-B	Brown	RS485-B	Blue
O	RS485-A	Pink	RS485-A	White
P	+Ain ⁴⁾	–	+Ain ⁴⁾	–
Q	–Ain ⁴⁾	–	–Ain ⁴⁾	–

1) PROFIBUS PA® only

2) Connect the grey +24 DC wire and the blue GND_24 V wire of the sensor separately to an external power supply.

3) Connect the brown +24 DC wire and the black GND_24 V wire of the sensor separately.

4) 4 to 20 mA signal for pressure compensation

Ordering information

Transmitter	Order no.	Description
M400/2H, 1-channel multi-parameter	30 025 514	4 to 20 mA (with HART®), Non-Ex version
M400/2XH, 1-channel multi-parameter	30 025 515	4 to 20 mA (with HART®), Ex version
M400G/2XH, 1-channel multi-parameter	30 025 516	4 to 20 mA (with HART®), for gas applications, Ex version
M400/2XH Cond Ind, 1-channel	30 256 307	4 to 20 mA (with HART®), for analog inductive conductivity sensors, Ex version
M400FF, 1-channel multi-parameter	30 026 616	FOUNDATION fieldbus®, Ex version
M400PA, 1-channel multi-parameter	30 026 617	PROFIBUS PA®, Ex version

M400/2H, M400/2XH, M400G/2XH parameter fit guide

	Analog	ISM	Analog	ISM
pH/ORP	•	•	•	•
pH/pNa	–	•	–	•
Conductivity 2-e	•	–	•	–
Conductivity 4-e	•	•	•	•
Amp. dissolved oxygen ppm/ppb/trace	•/•/•	•/•/•	•/•/•	•/•/•
Amp. oxygen gas	–	–	•	•
Optical dissolved oxygen ppm/ppb	–	•/•	–	•/•
Dissolved carbon dioxide (CO ₂ low)	–	•	–	•

M400FF, M400PA parameter fit guide

	Analog	ISM	Analog	ISM
pH/ORP	•	•	•	•
pH/pNa	–	•	–	•
Conductivity 2-e	•	–	•	–
Conductivity 4-e	•	•	•	•
Amp. dissolved oxygen ppm/ppb/trace	•/•/•	•/•/•	•/•/•	•/•/•
Amp. oxygen gas	•	•	•	•
Optical dissolved oxygen ppm/ppb	–	•/•	–	•/•
Ozone	–	–	–	•
Dissolved carbon dioxide (CO ₂ low)	–	•	–	•
Thermal conductivity (CO ₂ hi) (InPro 5500 i)	–	–	–	•

M400/2XH Cond Ind parameter fit guide

	Analog
Cond Ind (Inductive Conductivity) ¹⁾	•

1) InPro 7250 ST, InPro 7250 PFA. InPro 7250 HT

Accessories

Description	Order no.
Pipe Mount Kit for ½DIN models	52 500 212
Panel Mount Kit for ½DIN models	52 500 213
Wall Mount Kit for ½DIN models	30 300 482
Protective Hood	52 500 214

Sales and Service:

Australia

Mettler-Toledo Limited
220 Turner Street
Port Melbourne, VIC 3207
Australia
Phone +61 1300 659 761
e-mail info.mtaus@mt.com

Austria

Mettler-Toledo Ges.m.b.H.
Laxenburger Str. 252/2
AT-1230 Wien
Phone +43 1 607 4356
e-mail prozess@mt.com

Brazil

Mettler-Toledo Ind. e Com. Ltda.
Avenida Tamboré, 418
Tamboré
BR-06460-000 Barueri/SP
Phone +55 11 4166 7400
e-mail mtbr@mt.com

Canada

Mettler-Toledo Inc.
2915 Argenta Rd #6
CA-ON L5N 8G6 Mississauga
Phone +1 800 638 8537
e-mail ProlnsideSalesCA@mt.com

China

Mettler-Toledo International Trading
(Shanghai) Co. Ltd.
589 Gui Ping Road
Cao He Jing
CN-200233 Shanghai
Phone +86 21 64 85 04 35
e-mail ad@mt.com

Croatia

Mettler-Toledo d.o.o.
Mandlova 3
HR-10000 Zagreb
Phone +385 1 292 06 33
e-mail mt.zagreb@mt.com

Czech Republic

Mettler-Toledo s.r.o.
Trebohosticka 2283/2
CZ-100 00 Praha 10
Phone +420 2 72 123 150
e-mail sales.mtcz@mt.com

Denmark

Mettler-Toledo A/S
Naverland 8
DK-2600 Glostrup
Phone +45 43 27 08 00
e-mail info.mtdk@mt.com

France

Mettler-Toledo
Analyse Industrielle S.A.S.
30, Boulevard de Douaumont
FR-75017 Paris
Phone +33 1 47 37 06 00
e-mail mtpro-f@mt.com

Germany

Mettler-Toledo GmbH
Prozeßanalytik
Ockerweg 3
DE-35396 Gießen
Phone +49 641 507 444
e-mail prozess@mt.com

Great Britain

Mettler-Toledo LTD
64 Boston Road, Beaumont Leys
GB-Leicester LE4 1AW
Phone +44 116 235 7070
e-mail enquire.mtuk@mt.com

Hungary

Mettler-Toledo Kereskedelmi KFT
Teve u. 41
HU-1139 Budapest
Phone +36 1 288 40 40
e-mail mth@axelero.hu

India

Mettler-Toledo India Private Limited
Amar Hill, Saki Vihar Road
Powai
IN-400 072 Mumbai
Phone +91 22 2857 0808
e-mail sales.mtin@mt.com

Indonesia

PT. Mettler-Toledo Indonesia
GRHA PERSADA 3rd Floor
Jl. KH. Noer Ali No.3A,
Kayuringin Jaya
Kalimalang, Bekasi 17144, ID
Phone +62 21 294 53919
e-mail
mt-id.customersupport@mt.com

Italy

Mettler-Toledo S.p.A.
Via Vialba 42
IT-20026 Novate Milanese
Phone +39 02 333 321
e-mail
customercare.italia@mt.com

Japan

Mettler-Toledo K.K.
Process Division
6F Ikenohata Nissshoku Bldg.
2-9-7, Ikenohata
Taito-ku
JP-110-0008 Tokyo
Phone +81 3 5815 5606
e-mail helpdesk.ing.jp@mt.com

Malaysia

Mettler-Toledo (M) Sdn Bhd
Bangunan Electroscon Holding, U 1-01
Lot 8 Jalan Astaka U8/84
Seksyen U8, Bukit Jelutong
MY-40150 Shah Alam Selangor
Phone +60 3 78 44 58 88
e-mail
MT-MY.CustomerSupport@mt.com

Mexico

Mettler-Toledo S.A. de C.V.
Ejército Nacional #340
Polanco V Sección
C.P. 11560
MX-México D.F.
Phone +52 55 1946 0900
e-mail mt.mexico@mt.com

Norway

Mettler-Toledo AS
Ulvenveien 92B
NO-0581 Oslo Norway
Phone +47 22 30 44 90
e-mail info.mtn@mt.com

Poland

Mettler-Toledo (Poland) Sp.z.o.o.
ul. Poleczki 21
PL-02-822 Warszawa
Phone +48 22 545 06 80
e-mail polska@mt.com

Russia

Mettler-Toledo Vostok ZAO
Sretenskij Bulvar 6/1
Office 6
RU-101000 Moscow
Phone +7 495 621 56 66
e-mail inforus@mt.com

Singapore

Mettler-Toledo (S) Pte. Ltd.
Block 28
Ayer Rajah Crescent #05-01
SG-139959 Singapore
Phone +65 6890 00 11
e-mail
mt.sg.customersupport@mt.com

Slovakia

Mettler-Toledo s.r.o.
Hattalova 12/A
SK-831 03 Bratislava
Phone +421 2 4444 12 20-2
e-mail predaj@mt.com

Slovenia

Mettler-Toledo d.o.o.
Pot heroja Trtnika 26
SI-1261 Ljubljana-Dobrunje
Phone +386 1 530 80 50
e-mail keith.racman@mt.com

South Korea

Mettler-Toledo (Korea) Ltd.
1 & 4F, Yeil Building 21
Yangjaecheon-ro 19-gil
SeoCho-Gu
Seoul 06753 Korea
Phone +82 2 3498 3500
e-mail Sales_MTKR@mt.com

Spain

Mettler-Toledo S.A.E.
C/Miguel Hernández, 69-71
ES-08908 L'Hospitalet de Llobregat
(Barcelona)
Phone +34 902 32 00 23
e-mail mtemkt@mt.com

Sweden

Mettler-Toledo AB
Virkesvägen 10
Box 92161
SE-12008 Stockholm
Phone +46 8 702 50 00
e-mail sales.mts@mt.com

Switzerland

Mettler-Toledo (Schweiz) GmbH
Im Langacher, Postfach
CH-8606 Greifensee
Phone +41 44 944 47 60
e-mail ProSupport.ch@mt.com

Thailand

Mettler-Toledo (Thailand) Ltd.
272 Soi Soonvijai 4
Rama 9 Rd., Bangkok
Huay Kwang
TH-10320 Bangkok
Phone +66 2 723 03 00
e-mail
MT-TH.CustomerSupport@mt.com

Turkey

Mettler-Toledo Türkiye
Haluk Türksoy Sokak No: 6 Zemin ve 1.
Bodrum Kat 34662 Üsküdar-İstanbul, TR
Phone +90 216 400 20 20
e-mail sales.mtr@mt.com

USA

METTLER TOLEDO
Process Analytics
900 Middlesex Turnpike, Bld. 8
Billerica, MA 01821, USA
Phone +1 781 301 8800
Freephone +1 800 352 8763
e-mail mtprous@mt.com

Vietnam

Mettler-Toledo (Vietnam) LLC
29A Hoang Hoa Tham Street, Ward 6
Binh Thanh District
Ho Chi Minh City, Vietnam
Phone +84 8 35515924
e-mail
MT-VN.CustomerSupport@mt.com



Management System
certified according to
ISO 9001 / ISO 14001

Subject to technical changes.
02/2017. © Mettler-Toledo GmbH
Printed in Switzerland. 30 031 507

Mettler-Toledo GmbH, Process Analytics
Im Hackacker 15, CH - 8902 Urdorf, Switzerland
Phone + 41 44 729 62 11, Fax +41 44 729 66 36

www.mt.com/pro