


InLab[®] 742-ISM

Cond °C ISM



For portable applications in the lab and outdoor. Ideal for low conductivity range samples like pure water and ultra pure water.



Sensor type conductivity cell

Order number 51344116

Specifications

measuring range 0.001...500 $\mu\text{S}/\text{cm}$

zero point (pH) -

temperature range 0...100 °C

temperature probe NTC 30 k Ω

shaft material stainless steel

membrane glass -

membrane resistance -

reference system -

reference diaphragm -

reference electrolyte -

bridge electrolyte -

electrode head -

cable length 1.8 m

connector LTW (IP67)

cell constant (cm^{-1}) 0.105 (certified)

cell type 2 steel poles

indicating metal -

Storage

dry