


# InLab<sup>®</sup> Micro Pro-ISM

pH Ref °C ISM



Micro sensor for small sample volumes like test tubes, vials

Scale 1:1



**Sensor type** combined pH electrode

**Order number** 51344163

### Specifications

pH range 0...14

zero point (pH)  $7.00 \pm 0.25$

temperature range 0...100 °C

temperature probe NTC 30 k $\Omega$

shaft material glass

membrane glass U

membrane resistance < 300 M $\Omega$

reference system ARGENTHAL<sup>™</sup>

reference diaphragm ceramic

reference electrolyte 3 mol/L KCl

bridge electrolyte -

electrode head MultiPin<sup>™</sup>

cable length -

connector -

cell constant (cm<sup>-1</sup>) -

cell type -

indicating metal -

### Storage

InLab storage solution (30111142)

or 3 mol/L KCl (51350072)