

## Simple Data Management EasyDirect™ Balance Software



### Collect data automatically

Collect weighing data from up to 10 balances via Ethernet, USB-C or RS232. Minimize manual transcription errors, save time and store data securely.



### Efficient results handling

Easily review results and filter by date, instrument, user or sample. Visualize results in charts to assess target and tolerance ranges and perform statistics for efficient trend and production analysis.



### Easy reporting and export

Generate simple and clear reports of weighing results. Export data to PC in various formats (XML, CSV, XLSX or PDF) or print on a network printer.



### Simple instrument overview

By providing adjustment, test and service status information, EasyDirect Balance gives you a simple operational overview of all your connected instruments.



### Weighing Data-Handling Made Easy

Connect your balances to EasyDirect Balance PC software and improve the management and handling of your data. Installation and setup are simple and fast, and operation is intuitive.

Easy data management will help improve the efficiency of your weighing tasks. With EasyDirect Balance, you collect weighing data automatically, making it easier to create reports, visualize results in control charts, perform statistics and have a clear overview about the status of your balances.

Download your free 45-day trial version at:

► [www.mt.com/EasyDirectBalance](http://www.mt.com/EasyDirectBalance)

## Technical Description

Supported Instruments	MS-TS, ML-T <sup>1)</sup>	MX, MR, MA, LA, MS, ML, ME-T, ME, PL-E, JP, JS, JE, JET, JL-GE, XP, XPE, XS, XSE
<b>Connectivity</b>	Ethernet, RS232	RS232, USB-C
<b>Data Collection</b>		
<b>Balance Information</b>		
Basic balance information	•	•
Level, adjustment, test and service status	•	
Balance usage	•	
<b>Weighing Data</b>		
Supported applications	Weighing, counting, percent weighing, factor weighing	•
	Checkweighing, dynamic weighing, formulation, totaling, back weighing, differential weighing, density	•
Weight, unit, date/time, sample and task IDs	•	•
Target and tolerances	•	
User name	•	
Application-specific results and parameter	•	
<b>Maximum Number of Instruments</b>	Up to 10 balances	
<b>Report Formats</b>	XML, CSV, XLSX, PDF	
<b>Database</b>	Microsoft SQL Server 2012, Express Edition	
<b>PC Requirements</b>	OS: Windows 10 and 11	
<b>Protection</b>	Password protection of settings and results	
<b>Languages</b>	English, German, French, Spanish, Italian, Russian, Polish, Czech, Hungarian, Dutch, Portuguese, Turkish, Japanese, Chinese, Korean	
<b>Material Numbers</b>	30539323 EasyDirect Balance (License for up to 3 instruments) 30540473 EasyDirect Balance (License for up to 10 instruments)	

<sup>1)</sup> Full functionality only via Ethernet



**METTLER TOLEDO Group**  
 Laboratory Weighing  
 Local contact: [www.mt.com/contacts](http://www.mt.com/contacts)

[www.mt.com/EasyDirectBalance](http://www.mt.com/EasyDirectBalance)

For more information

Subject to technical changes  
 © 03/2024 METTLER TOLEDO. All rights reserved  
 30553484B  
 Group MarCom RITM1159092 MB/NM